

Rabbit IgG anti-Sheep IgG (Fc)-Alk. Phos., MinX Hu

General Information

Catalog Number:	SEC-183559	Physical State:	Liquid (sterile filtered)
Quantity:	0,5 mg	Reconstitution Volume:	
		Reconstitution Buffer:	
		Shipping Conditions:	Wet Ice
		Product Expiration	Expiration date is one (1) year from date of opening.

Antibody Host / Format

Host Species:	Rabbit	Clonality / Isotype / Clone:	polyclonal Ig
Antibody Format:	IgG whole molecule	Concentration:	1.0 mg/mL by UV absorbance at 280 nm
		Conjugation:	Alkaline Phosphatase
		Maximum Excitation / Emission:	/

Specificity

Target Species:	Sheep	Immunogen:	Sheep IgG F(c) fragment
Antibody Specificity:	IgG (Fc Fragment)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Sheep IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Alkaline Phosphatase (calf intestine), anti-Rabbit Serum, Sheep IgG, Sheep IgG F(c) and Sheep Serum. No reaction was observed against Sheep IgG F(ab) or Human Serum Proteins.. Minimal Cross Reactivity: Human	

Formulation, Transport and Storage

Storage Buffer:	0.05 M Tris Chloride, 0.15M Sodium Chloride, 0.001M Magnesium Chloride, 0.0001M Zinc Chloride, 50% (v/v) Glycerol; pH 8.0	Preservatives:	0.01% (w/v) Sodium Azide
Stabilizer:	10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free	Store vial at 4° C before opening. DO NOT FREEZE. This product is stable at 4° C as an undiluted liquid. Dilute only prior to immediate use. Freezing alkaline phosphatase conjugates will result in a substantial loss of enzymatic activity. Storage Temperature: 4-8°C	

Application Recommendation

Manufacturer

BIOZOL
Diagnostica Vertrieb GmbH
Oehleckerring 11-13
22419 Hamburg

Phone +49 (89) 3799 666 6
Fax +49 (89) 3799 666 99
E-mail info@biozol.de

Find Secondary Antibodies

www.biozol.de



Application: ELISA, Immunohistochemistry (IHC), Immunohistochemistry (Paraffin-embedded Sections), Immunohistochemistry (frozen sections), Western Blot

Dilution: ELISA 1:2,000 - 1:10,000, Immunohistochemistry (IHC) 1:200 - 1:1,000, Western Blot (WB) 1:500 - 1:2,500

Application Note:

This product has been assayed against 1.0 ug of Sheep IgG in a standard capture ELISA using pNPP p-nitrophenyl phosphate code # NPP-10 as a substrate for 30 minutes at room temperature. A working dilution of 1:3,000 to 1:14,000 of the reconstitution concentration is suggested for this product.

Background Information

Anti-Sheep IgG F(c) generated in rabbit is a proteolytic fragment of immunoglobulin G (IgG) obtained by limited digestion with the enzyme papain under controlled conditions of temperature, time and pH. Receptors bind the Fc portion of sheep IgG and often this fragment is removed from immunoglobulins to minimize receptor binding and lower background reactivity.

Safety Information

This reagent contains sodium azide (a poisonous and hazardous substance) as preservative. Although this concentration is not regarded as dangerous to health, appropriate measures should be taken to avoid skin and eye contact, inhalation, and ingestion. Standard Laboratory Practices should be followed when handling this material.

Disclaimer

This product is warranted to perform in conformance with product specifications and to be free from defects in material and workmanship. Products are supplied for research use only. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated, this warranty is limited to one year from date of shipment when the product is subjected to normal, proper and intended usage. Buyer's exclusive remedy for non-conforming products during the warranty period is limited to replacement of or refund for the non-conforming product(s). There is no obligation to replace products as the result of accident or disaster or inappropriate use of the products in a manner for which they were not designed, or improper storage and handling of the products.

Manufacturer

BIOZOL
Diagnostica Vertrieb GmbH
Oehleckerring 11-13
22419 Hamburg

Phone +49 (89) 3799 666 6
Fax +49 (89) 3799 666 99
E-mail info@biozol.de

Find Secondary Antibodies

www.biozol.de

